

ENERGY STATES OF ATOMS – THE FRANCK-HERTZ EXPERIMENT

Atoms in excited states emit radiation at discrete frequencies. The frequency of the radiation, ν , is related to the change in the atomic energy levels by $\Delta E = h\nu$. The fact that only certain discrete frequencies are emitted is evidence of the quantization of atomic energy levels. In the present experiment, this quantization is examined by bombarding mercury vapour with electrons of known, variable energy. This experiment is a repeat of one originally performed by Franck and Hertz in 1914. They set out to show that: i) it is possible to excite atoms by low-energy electron bombardment; ii) the excitation energy transferred from the electrons to the atoms always has discrete values; and iii) the values obtained for the energy levels are in agreement with spectroscopic results.

Theory:

When particles travel through matter there are two possible types of collisions, elastic and inelastic, between the incident particles and those forming the substance. In an elastic collision kinetic energy is conserved and the initial energy is divided between the particles according to the relation of their masses.

If an electron (small mass) collides elastically with an atom (large mass), from conservation of momentum and kinetic energy it is seen that the electron's direction of motion is changed but that the amount of kinetic energy lost by the electron and imparted to the atom is very small.



If an electron collides inelastically with an atom, the energy lost by the electron is not converted into kinetic energy of the atom, but rather is the cause for excitation, ionization, radiation, or particle production. Quantum theory postulates that an atom can only be excited in discrete energy steps, i.e. there is a small but finite energy required to bring an atom from its ground state into an excited state. If, in a collision, the energy available to be transferred to the atom is less than the difference between the ground and first excited states, the collision remains elastic. However, an electron having a kinetic energy larger than the energy difference between the ground state and the first excited state of an atom can excite this atom by means of an inelastic collision.

There are other ways to excite atoms. Illumination with light of the appropriate wavelength (given by $\lambda = hc/\Delta E$, where ΔE = energy difference between ground and an excited state) would cause a transition. The inconvenience of using light to excite an atom is that this light must have the exact wavelength corresponding to the energy difference between the ground state and the state to be excited. In using the energy loss due to an inelastic collision as the energy transport mechanism, this drawback does not occur since the kinetic energy of the electron is not quantized either before or after the collision.

Apparatus:

There are two sets of equipment. One set enables examination of the collisions between electrons and mercury atoms and the other set of equipment enables examination of the collisions between electrons and neon atoms.


Franck-Hertz Experiment in Mercury





The Franck-Hertz tube is an evacuated glass cylinder containing a drop of mercury. The tube contains three electrodes: indirectly-heated cathode, grid-form anode, and collector electrode. The electrodes are arranged in plane-parallel manner.

The entire tube is housed in an oven enclosure. When the tube is heated in the oven, some of the mercury vaporises. After about 20 minutes the temperature reaches equilibrium around 200 °C. The Hg vapour pressure is now high enough that the mean free path of the electrons emitted by the cathode is smaller than the distance between cathode and anode.




The retarding voltage, V_r , is provided by a standard AA 1.5 V battery.

The accelerating voltage, V_a , is provided by an ELWE Franck-Hertz Operating Unit.



An Electro Industries DC Power Supply is used to supply the current to heat the cathode so that electrons are emitted.



A Keithley 610B Electrometer is connected to the counter electrode to measure the electron current in the tube. The output terminals at the back of the electrometer, which provide a voltage that is proportional to the detected current, are connected to a digital voltmeter.

Data is collected manually by reading the accelerating voltage that corresponds to each of the current peaks, as observed by watching the electrometer meter.

Franck-Hertz Experiment in Neon


The equipment used to study the electron-atom interactions in neon is similar to that used for the Franck-Hertz in Mercury experiment, except that an oven is not required and the Franck-Hertz Operating Unit is used to provide all the necessary voltages. The operating unit also contains a high-sensitivity DC amplifier for measuring the electron current at the counter electrode.



Data is collected using a PASCO 550 computer interface. One input of the interface is used to collect the accelerating voltage data from the $U_B/10$ output of the F-H Operating Unit and the other input of the interface is used to collect electrometer current data.



The operating unit provides outputs for displaying the Franck-Hertz signal on a dual-channel oscilloscope.




Characteristics of the Electron-Atom Interaction

The emitted electrons gain energy due to the accelerating voltage V_a applied between cathode and anode. After passing through the perforated anode, the electrons encounter a retarding, decelerating voltage V_r . Since V_r is about 1.5 V, electrons having an energy of less than 1.5 eV are repelled and do not contribute to the current flowing from the counter-electrode to the electrometer. The current, as measured at the electrometer, indicates how many electrons passed the anode and reached the counter-electrode. If V_a is smaller than the energy difference between the ground state and first excited state of the atoms in the tube, the collisions between electrons and the atoms remain elastic, so the electrons change their direction of flight but lose very little kinetic energy due to the impact. The electron current rises with increasing V_a in accordance with normal tube characteristics, proportional to $V_a^{3/2}$. If V_a is further increased, the energy of the electrons becomes sufficient to raise the atoms into the first excited state during a collision. That is, the collisions become inelastic. Since the impacting electron loses this excitation energy, the remaining kinetic energy will not be sufficient for the electron to overcome the retarding potential V_r and consequently this electron is not detected. Thus the current will drop directly after V_a exceeds the excitation energy, sometimes called the critical or resonance energy.

If the accelerating voltage is raised further, the electron current increases since the electrons can gain enough energy after the inelastic collision to reach the counter-electrode. If V_a exceeds twice the excitation energy, electrons will undergo a second inelastic collision and consequently the current will drop once more. This pattern will repeat as V_a is increased.

The following figure shows the expected plot of electrometer current I versus accelerating voltage V_a .



Mercury Source – Detailed Procedure:

1. Plug in the oven, if this has not already been done, and allow approximately 20 minutes for the temperature to reach equilibrium before taking measurements.
2. Check that the Keithley 610B Electrometer is set as follows:

Range control	3×10^{-10} A
FEEDBACK switch	FAST
ZERO CHECK switch	LOCK position (push and turn)

Turn on the electrometer by turning the Meter switch to the “-” position and allow five to ten minutes for warm-up. Unlock the ZERO CHECK switch.

3. The cathode current is supplied by the Electro Industries Regulated DC Power Supply. Check that the VOLTAGE RANGE switch is set to 0-15 V, that the current control is at minimum (fully counterclockwise), and that the voltage control is set at midrange (indicator on knob pointing vertically upward). Turn on the digital multimeter connected to this supply by setting it to the 20m/10A DCA position. Turn on the power supply. **Slowly** increase the current until the ammeter reads about 0.36 A.
4. The accelerating voltage is provided by the Franck-Hertz Operating Unit. Set the ACCELERATION control to 0 and put the toggle switch in the Man. (manual) position. Connect the A (anode) terminal of the operating unit to the anode connection on the tube cabinet. Connect the K (cathode) terminal of the operating unit to the rightmost cathode

terminal on the tube cabinet. Turn on the multimeter that is connected across the anode-cathode terminals on the tube cabinet and set the meter to display DC voltage.

5. Turn on the multimeter that is connected to the back of the electrometer. Set the meter to the 200 DCV position. This multimeter measures a voltage that is proportional to the number of electrons reaching the counter electrode.
6. The system is now ready for use.

By turning the ACCELERATION control on the Franck-Hertz Operating Unit, slowly and smoothly increase the accelerating voltage from 0 to about 25 volts, while watching the tube current on the meter of the Keithley electrometer.

Record the values of accelerating voltage corresponding to the peaks in tube current.

When near a peak, you will want to vary the voltage back and forth over a small range in order to accurately determine the voltage that corresponds to the current peak.

When finished recording your data, turn off the equipment as follows:

Operating Unit:	Set voltage to 0 Turn off unit
Cathode Current Supply:	Reduce cathode current to 0 Turn off supply
Electrometer:	ZERO CHECK to LOCK position (push and turn) POWER OFF
Oven:	ask the instructor to unplug the cord from the power bar
Multimeters:	switch to OFF

Mercury Source – Analysis:

From your recorded Accelerating Voltage values corresponding to each of the current peaks, determine an average value for the peak separations, and hence for the energy difference between the ground state and the first excited state of mercury. Based on your experimental value for the energy of the first excited state of mercury, calculate the wavelength of light that would be emitted when mercury de-excites from the first excited state to the ground state. Compare with the accepted value of 253.7 nm.

What do you notice about the location of the first peak compared to the average peak separation? Attempt to explain this discrepancy.

Why does the current decrease gradually, rather than suddenly, at the critical energies?

Do you expect to see effects due to excitation of the second and higher excited states of mercury? Explain your answer.

Neon Source – Detailed Procedure:

1. If necessary, unplug the leads connecting the A and K terminals of the Franck-Hertz Operating Unit to the F-H Mercury tube.
2. Connect the leads from the F-H Neon tube to the appropriate terminals of the operating unit.
3. If necessary, disconnect the PASCO 550 CH A voltage sensor leads from the electrometer output voltmeter and connect them to the FH signal terminals of the operating unit.
4. Connect the coax signal cable from the neon tube to the FH Signal input of the operating unit.
5. Connect CH1 of the oscilloscope to the $U_B/10$ output terminals of the operating unit and connect CH2 of the oscilloscope to the FH signal output terminals.
6. Set the operating unit controls and the oscilloscope as follows:

Reverse Bias:

- ◆ just over 4 (just past indicator mark between 2 and 6)


Acceleration Voltage:

- ◆ ~65 V
- ◆ switch set to Ramp 50 Hz

Heater Voltage:


- ◆ ~8 V
- ◆ tube behaviour and output trace are strongly dependent on heater voltage
- ◆ ***Once the trace appears on the oscilloscope, you may need to slightly reduce the heater voltage***

Amplitude:


- ◆ set about halfway between \longrightarrow and  (about 25% of full scale)

Oscilloscope Settings (x-y mode):

- ◆ X: 1 V/cm Y: 0.5 V/cm
 - ◆ Note that the horizontal output of the control unit is one-tenth of the accelerating voltage.
7. Turn on the F-H Operating Unit and the oscilloscope and allow approximately 10 minutes for the equipment to stabilize.
 8. Adjust the oscilloscope scales as necessary to obtain a stable trace showing the neon excitation peaks:



9. From the oscilloscope display, measure the accelerating voltages corresponding to as many excitation peaks as are visible. **Note that the voltage along the x -axis of the oscilloscope display is one-tenth of the actual accelerating voltage that is being applied to the F-H tube.**
10. Decrease the accelerating voltage to zero. Slowly increase the accelerating voltage while carefully observing the tube through the viewing port. **You will need to dim the room lights.** The following image shows the F-H Neon tube in operation.




11. Record your observations in as much detail as possible. In particular, ensure you have enough information to answer the following questions. What colour is the light produced in the tube? Where is the light first produced? Does the location of the light depend on accelerating voltage? If so, how? Do you observe bands of light in the tube?
12. Data can also be acquired, analysed, and saved using the PASCO 550 interface.
13. On the Franck-Hertz Operating Unit, set the ACCELERATION control to 0 and set the toggle switch to Man.
14. Channel A of the PASCO 550 is connected to a CI-6503 voltage sensor. The leads from this voltage sensor are connected to the FH Signal and ground outputs of the Franck-Hertz Operating Unit.

15. Channel B of the PASCO 550 is also connected to a CI-6503 voltage sensor. The leads from this voltage sensor are connected to the $U_B/10$ and ground outputs of the Franck-Hertz Operating Unit.
16. Turn on the PASCO 550 interface.
17. Double-click the PASCO Capstone icon on the desktop of the computer.
 - Click the Table and Graph template (the first option)
 - At the top of the Tools menu on the left of the window, click Hardware Setup
 - On the graphic that opens, click CH A, scroll down the list of sensors, and click Voltage Sensor
 - Click CH B, scroll down the list of sensors, and click Voltage Sensor
 - Close the Hardware Setup window by clicking on the Hardware Setup icon
 - At the bottom middle of the window, set the Voltage Sensor sample rate to 20.00 Hz
 - Ensure that both voltage sensors are set to a sample rate of 20.00 Hz
 - In the data table, click <Select Measurement> in the left-hand column and select Voltage Ch A
 - Do the same for the right-hand column, but select Voltage Ch B under Voltage Sensor
 - On the graph, set the y-axis to display the Ch A Voltage and the x-axis to display the Ch B Voltage by clicking on the <Select Measurement> buttons.
 - Label the graph axes as F-H Current for the vertical axis and Accelerating Voltage for the horizontal axis. This is done by clicking near the appropriate axis and selecting the gear icon from the toolbar above the graph. Click Axis Label, Show Custom Label, and Edit the Custom Label appropriately.
18. The equipment is now ready for use.
19. Click the Record button in the Capstone software and slowly and steadily increase the accelerating voltage by turning the ACCELERATION control on the operating unit.
20. When the accelerating voltage reaches about 70 V click the Stop button in Capstone. It should take 45 to 60 seconds to reach 70 V.
21. As the interface records the data, it will appear in both the data table and graph within Capstone.
22. The data can either be analysed within Capstone or saved to a file for later analysis using a spreadsheet program such as Microsoft Excel.
23. To analyse the data within Capstone, note that the graph can be easily manipulated:
 - clicking and dragging the plot area near either of the axes allows adjustment of the starting values on the graph axes;
 - clicking and dragging the numbers on each of the axes allows the axes scales to be adjusted independently.
24. To save the data as a tab-delimited text file, click File, Export Data.... Choose an appropriate location for the data file, give it a descriptive filename, and click Save. The data file can be opened in Excel for graphing and analysis at a later time.

To delete data, use the Delete Last Run or Delete All Runs button at the bottom right of the window (in the Controls toolbar).

Graph of typical data:



25. When finished recording your data, turn off the oscilloscope, reduce the accelerating voltage to 0, and turn off the Franck-Hertz Operating Unit.
26. Export your data, then close the Capstone program.
27. Turn off the PASCO 550.

Neon Source – Analysis:

Determine an average value for the peak separations, and hence for the energy difference between the ground state and the excited state of neon. Remember that the output of the accelerating voltage terminals of the F-H Operating Unit is one-tenth of the actual accelerating voltage that was applied to the tube.

Calculate the wavelength of light that would be emitted when neon de-excites from the excited state to the ground state.

Is this wavelength in the appropriate part of the visible spectrum to account for the colour of light that is observed to be emitted by the excited neon?

If not, note that there is a band of neon states at ~ 16.7 eV. Calculate the wavelength of light that would be emitted for transitions between the excited state that you measured and the states at 16.7 eV. Does this wavelength correspond to the colour of light emitted by the neon?

References:

- Melissinos, Experiments in Modern Physics, QC 33
- Harnwell, Experimental Atomic Physics, QC 173